



**HAIGH**

Introducing

# Panaway<sup>®</sup> M1

The first movable bedpan disposer

## **INFECTION PREVENTION: MOVING FROM REUSABLE TO DISPOSABLE COMMODE POTS, BEDPANS AND URINE BOTTLES**

It is estimated 95% of NHS healthcare sites have already made the switch to disposable (single use) systems. More and more elderly care companies are enjoying the benefits every year:

- Globally recognised as part of best practice to preventing infections and reducing cross-contamination relating to human waste
- Quick and easy to use – giving staff more time for high-value activities
- Hygienic and environmentally friendly – pulp products are clean, biodegradable and made from sustainable materials
- Economical and energy efficient – using just cold water and minimal energy per cycle

If you are looking to move from washing and sterilising your commode pots to using disposable items, we can help.

### **An evolutionary step forward**

Haigh's Panaway<sup>®</sup> M1 is your perfect solution: the first movable bedpan macerator of its kind, leading the way in performance, versatility, and ease of use – with hygiene, economy and efficiency at the heart of its design. Available on a short-term trial basis and as part of a longer-term migration.



## Reliability you can trust

We know transitions can be tricky and that's why we offer a stress-free trial of our Panaway® M1, enabling you to:

- Test the macerator and familiarise staff with its functionality
- Gather performance data to inform asset maintenance and management decisions
- Enjoy low-risk rental agreements and economical pulp item supply

## Convenience without compromise

Panaway® M1 gives you the accessibility and benefit of a modern sluice room where – and when – you need it most, without compromising on safety.

- Compact and versatile, it can fit – and be wheeled – into the most convenient location for your site
- Adding measurable value as an emergency stand-in when a bedpan machine (washer or macerator) is down and awaiting maintenance
- A temporary solution whilst you migrate your site to disposable

## Simpler, safer infection prevention

Disposable pulp items and the contents\* can be placed into your conveniently located Panaway® M1.

- Easy to set up, Panaway® M1 is supplied with all fixtures and fittings required to use in a wide range of locations
- Fast, efficient and quiet, the Panaway® M1 disposes of the pulp product and contents\* in less than a minute
- Easy to clean with antimicrobial plastic surfaces helping you maintain hygiene standards

*\* the macerator is designed to process one standard pulp item per cycle containing liquids, human waste, tissue and macerator- friendly wipes*



## Easy to use, economical and energy efficient

Discreet and easy to use, with clear visual instructions, the Panaway® M1 aims to make the life of the user as straightforward as possible.

- Panaway® M1 uses less than a quarter of the energy and water of washer steriliser systems
- No hot water or heating is not just good value, it is good for the environment
- Paper pulp is eco-friendly, hygienic and highly cost-effective
- No annual / bi-annual heating element temperature test, mean low operating costs

Even if you do not wish to switch to disposables, a Panaway® M1 is a great temporary backup for when your washer is down for maintenance.



*The Panaway trial scheme was exactly what we needed to test the transition to disposable pulp items. I'm thrilled – it's facilitated our future!"*

*Manager, Private Care Home*



## Benefits of Panaway® M1

### FOR PATIENTS AND CARERS

- Reduces risk of infection and cross-contamination
- Supports discreet, quiet and fast waste disposal
- Requires minimal training as straightforward to use
- Frees up time for high value activities
- A great backup plan for when the washer is down

### FOR MANAGERS / OWNERS

- Proven reduction risk of infection and cross-contamination compared to washing
- Cost effective, energy efficient and environmentally friendly – no heating or hot water requirements
- Reliable, flexible, compact – trial layouts with low commitment
- Easy to set up and relocate (needs only a cold-water tap, 240v power supply and drainage)
- Low-risk trial and migration process



### Find out more

Contact your usual Haigh representative or **Michelle Marriott** on **+44(0)7500 626463** to discuss your specific situation and requirements.

For information about renting, purchasing or a trial of a Panaway® M1 email **info@haighmed.com**

**[www.haighmed.com/pm1](http://www.haighmed.com/pm1)**

**HAIGH**

