

# HAIGH



## Introducing Panaway® M1 The first movable bedpan disposer

### AN INFECTION PREVENTION STRATEGY

#### What do you do when your bedpan macerator breaks down?

Patient waste disposal management is an essential part of your healthcare facility's infection prevention and control strategy. A bedpan macerator breakdown is an immediate risk to those under your care.

Without an alternative method of disposal, hazardous waste products build up quickly, increasing infection risks for patients and staff. As a healthcare provider, the responsibility lies with you to provide a solution and fast.

#### An evolutionary step forward

With all the benefits and the proven reliability our customers have come to trust and expect from our products, we bring you Panaway® M1: the first movable bedpan macerator of its kind, leading the way in performance, versatility, and ease of use – with hygiene, economy and efficiency at the heart of its design.



## Convenience without compromise

Panaway® M1 gives you the accessibility and benefit of a sluice room where – and when – you need it most, without compromising on safety.

- Compact and versatile, it can fit – and be wheeled – into the most convenient location for your site
- Adding measurable value as an emergency stand-in when a macerator is down and awaiting maintenance
- A quick fix solution when awaiting installation of a fixed (Haigh) appliance

## Simpler, safer infection prevention

Disposable pulp products and their contents\* can simply be placed into your conveniently located Panaway® M1.

- Easy to set up, Panaway® M1 is supplied with all fixtures and fittings required to use in a wide range of locations
- Fast, efficient and quiet, the Panaway® M1 disposes of the pulp container and contents\* in less than a minute
- Easy to clean, with bacteria suppressing plastic and ergonomic design helps you maintain hygiene standards

## Easy to use, economical and energy efficient

Discreet and easy to use, with clear visual instructions, the Panaway® M1 aims to make the life of the user as straightforward as possible.

- Reduces stress associated with emergency breakdowns
- Panaway® M1 uses 35% less water than other macerators (and significantly less energy and water than washing systems)
- Paper pulp is eco-friendly, hygienic and highly cost-effective
- No annual / bi-annual heating element temperature test or hot water mean low operating costs

*\* the contents can include liquids, human waste, tissue paper and macerator-friendly wipes*



“

*“This is a game changer. Machine breakdowns and maintenance can be better managed and staff can focus on what they do best – caring for our patients.”*

*Senior Link Nurse, NHS Hospital*



## Benefits of Panaway® M1

### FOR PATIENTS AND CARERS

- Reduces risk of infection and cross-contamination
- Supports discreet, quiet and fast waste disposal
- Requires minimal training as straightforward to use
- Helps prevent overloaded sluice-rooms
- Reduces stress and dirty sluice rooms associated with breakdowns

### FOR MANAGERS / OWNERS

- Turns emergency callouts into less expensive, standard callouts – cost effective, energy efficient and environmentally friendly
- Reliable, flexible, compact – deploy where there is a need
- Easy to set up and relocate (needs only a cold-water tap, 240v power socket and drainage)
- No expensive maintenance or hot water required
- Saves time, freeing up front line staff for patient-focused work

---

### Find out more

Contact your usual Haigh representative or **Michelle Marriott** on **+44(0)7500 626463** to discuss your situation and requirements.

For information about renting or purchasing a Panaway® M1, email **info@haigh.co.uk** or visit **www.haighmed.com/pm1** to see our frequently asked questions.

**HAIGH**

